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- David A. Warrell, Timothy M. Cox, and John D. Firth Oxford Textbook of Medicine (5 ed.)
- Chetana Vaishnavi infections of the Gastrointestinal System
- Biosafety in Microbiological and Biomedical Laboratories, 5th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (2007).
- Henretig F.M. MD, King C. MD, Textbook of Pediatric Procedures, Chapter 123 – Obtaining Biologic Specimens Williams and Williams (April 1997).
- The Clinical Virology Laboratory, Department of Laboratory Medicine at Yale: http://info.med.yale.edu/labmed/virology/booklet.html.

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www.pbf.com.tw

www.vstriptech.com

Vstrip Adenovirus Rapid Test for the Detection of Adenovirus Antigen in Human Stool Specimen

INTENDED USE

Vstrip Adenovirus Rapid Test is an *in vitro* immunochromatographic assay for the qualitative detection of Adenovirus antigen in human stool specimen. The test is intended to aid in the early diagnosis of Adenovirus infection. This test is for professional use only.

SUMMARY AND EXPLANATION

Adenovirus is one of the most common cause of viral gastroenteritis in children. Depending on the serotype, this virus may also cause respiratory diseases, conjunctivitis and cystitis, etc. Adenovirus is non-enveloped, double-stranded DNA virus, 70-80 nm diameter, with icosahedral symmetry. At least 53 serotypes of Adenoviruses have been described in humans, all sharing a common hexon antigen. Serotypes 40 and 41 have been clearly associated with gastroenteritis disorders, whose main symptom is diarrhea, and may last for 5-9 days associated with temperature and vomits.^{1,2}

PRINCIPLE OF THE TEST

Vstrip Adenovirus Rapid Test is based on an immunochromatographic technology with colored microspheres. During testing, the sample is allowed to react with the colored conjugate (anti- Adenovirus monoclonal antibodies microspheres) which was pre-dried on the test strip. The mixture moves upward on the membrane by capillary action. The membrane is pre-coated with monoclonal antibodies, on the test band region, against viral antigens. If the sample extract contains Adenovirus antigens, these will form an antigen-antibody complex with the colored particles and a visible RED line will appear as the Test band (T). If Adenovirus antigen is not present, or is present at very low level in the specimen, only the BLUE line will be visible as the Control band (C). The BLUE line (C band) should always appear after a sample is applied, indicating a valid result. The device can accurately indicate the presence of Adenovirus in 10 minutes. If the control line does not develop, the test is invalid.

MATERIALS AND REAGENTS PROVIDED

1. Test dipsticks (20): Each strip contains Adenovirus specific monoclonal antibodies.
2. Stool collection tubes with diluent (20): The buffer contains detergent and 0.02% sodium azide.
3. Positive control reagent (1)
4. Package Insert (1)

MATERIALS NOT PROVIDED

1. Specimen collection container
2. Timer
3. Disposable gloves

WARNINGS AND PRECAUTIONS:

1. For *in vitro* diagnostic use.
2. Read the entire procedure carefully prior to testing.
3. Do not open the test pouches until ready to perform the assay.
4. Use of Nitrile or Latex gloves is recommended when handling patient samples.³
5. Do not use stools that have dried out.
6. Patient specimens should be mixed in the stool collection tubes with diluent gently and thoroughly before use.
7. Disregard test results beyond specified time (10 min).
8. Do not insert the test strip directly into patient stool specimens.
9. Do not reuse kit components or test devices.
10. Do not use the test if the foil pouch is damaged.
11. Do not mouth pipette samples or reagents.
12. Patient specimens may contain infectious agents and should be handled and disposed of as potential biohazards.
13. Do not interchange or mix different lots of Vstrip Adenovirus Rapid Test.
14. Do not use the kit components beyond the expiration date, printed on the outside of the box.
15. Do not insert the test strip directly into the mouth, eye or nasal.
16. Diluent buffer contains sodium azide, which may be a skin irritant. Avoid skin contact with reagents.
17. Dispose all used materials in the appropriate container. Treat as potentially biohazard.³
18. To obtain accurate results, you must follow the Package Insert.
19. Positive control reagent should be added 4-5 drops in the test tube and mixed gently before use.
20. Use appropriate precautions in the collection, handling, storage, and disposal of patient samples and used kit contents.³
21. Inadequate or inappropriate specimen collection, storage, and transport may yield false or negative test results.

22. Seek specific training or guidance if you are not experienced with specimen collection and handling procedures.^{4,5}
23. Specimens containing high levels of blood may fail to flow in the Vstrip Adenovirus Rapid Test device, resulting in an invalid/false positive test result.
24. Do not smoke, eat or drink in areas where specimens or kit components are handled.

STORAGE INSTRUCTION

The expiration date is indicated on the package label. Store kit reagents at 15-30°C, away from direct sunlight. The test must be kept in the sealed pouch until use. Do not freeze or overheat the test kit or kit reagents.

SPECIMEN COLLECTION AND HANDLING

Proper specimen collection, storage, and transport are critical to the performance of this test. For the best results, stool specimens should be collected after onset of symptoms. Do not leave specimens at room temperature for prolonged periods. Stool specimens should be collected in a clean container that do not contain media, preservatives, animal serum, or detergents as any of these additives may interfere with the Adenovirus antigen test. Specimens may be stored at 2-8°C for 2 days without interfering with the test performance. For long-term storage of specimens, -20°C or colder is recommended. Repeated freezing and thawing of specimens is not recommended and may cause erroneous results.

SPECIMEN PREPARATION

Stool specimen:

1. Unscrew the sample diluent vial cap, use the spoon attached to the cap to transfer small piece of stool (2-5 mm in diameter; approximately 30-50 mg or 30-50 ul) into the stool collection tube with diluent buffer.
2. Close the vial and tighten securely. Mix thoroughly by inversion, shaking or vortex in order to assure good sample dispersion.
3. Label the sample diluent vial with a patient or other control information.

ASSAY PROCEDURE

All clinical specimens must be at room temperature before beginning the assay.

Expiration date:

Check the expiration date on each individual test package or outer box before use. Do not use any test beyond the expiration date on the label.

1. Remove a test strip from a foil pouch.
2. Unscrew the cap of the stool collection tube.
3. Place the test strip into the stool collection tube with the arrows on the test strip pointing down. Do not handle or move the test strip until the test is complete and ready for interpretation.
4. Read the result at 10 minutes. Some strong positive results may appear sooner. Do not read result after 10 minutes.

INTERPRETATION OF RESULTS

Positive result:


A distinct RED colored line appears in the test region (T), in addition to a BLUE line appears in the control region (C). A positive result indicates the presence of Adenovirus antigen.

Negative result:

No RED line appears in the test region (T). A distinct BLUE line shows on the control line region (C). A negative result indicates that Adenovirus antigen is absent or below the level of detection.

Invalid:

No visually detectable the blue colored line in the control region, with or without a visually detectable red colored line on the test region. An excess of sample, incorrect procedure is mostly the main reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

QUALITY CONTROL

The Positive Controls should be assayed once upon receipt of the kit. Internal procedural controls are included in the test. A blue line appearing in the control region is an internal control. It confirms sufficient specimen volume and correct procedural technique. The positive controls are intended to monitor for reagent failure, but will not ensure precision at the analytical assay cutoff.

QUALITY CONTROL TESTING PROCEDURE

Positive Control:

1. Unscrew the positive control reagent and add 4-5 drops into the stool collection tube with diluent buffer.
2. Close the cap and mix the solution by inversion or shaking at least one time.
3. Remove the test strip from its foil pouch before using it.
4. The test strip put into the tube (arrows down).
5. Read the result at 10 minutes. Do not read result after 10 minutes.

LIMITATIONS OF THE PROCEDURE

1. A negative test result may occur if the level of antigen in a sample is below the detection limit of the test. If the symptoms or situation still present, a further determination should be carried out with other technique.
2. After one week of infection, the number of viruses in stool is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
3. A negative test result does not rule out other potential viral infections
4. An excess of sample could cause invalid result or may cause brown bands appearance. Dilute the sample with the diluent buffer and repeat the test.
5. Test Results must be evaluated in conjunction with other clinical data available to the physician.
6. A positive test result does not rule out co-infections with other pathogens or viruses.
7. A positive test result does not identify specific Adenovirus subtypes.

EXPECTED VALUES

Vstrip Adenovirus Rapid Test detects the presence of Adenovirus antigen in stool. Expected values for a given population should be determined for each laboratory. The prevalence of Adenovirus infection will vary based on many factors such as age, geographic location, method of sample collection, sample handling and transportation, and the general health environment of the patient population under study. Healthy individuals tested should be negative for Adenovirus. Some infected individuals may show symptoms or only minor symptoms, and these patients may test negative.

PERFORMANCE CHARACTERISTICS

The evaluation was conducted comparing the results obtained using Vstrip Adenovirus Rapid Test to another commercially available Adenovirus Rapid Test.

■ Sensitivity

The detection of Adenovirus showed >99% of concordance in sensitivity.

■ Specificity

The detection of Adenovirus showed > 99% of concordance in specificity.

■ Reproducibility

Three different lots of the Vstrip Adenovirus Rapid Test was tested in multiple replicates using negative, low positive, medium positive and high positive samples. Reproducibility was tested within each lot, between lots, and on different days. In each case, all tests yielded >99% reproducibility.

■ Limit of Detection

To determine the limits of detection of the Vstrip Adenovirus Rapid Test. The study used human Adenovirus as following table. Results are summarized below:

Adenovirus	Type	Limit of detection (LOD)
Adenovirus type	7	2.81 x 10 ⁴ (TCID)
Adenovirus type	41	1.25 x 10 ² (TCID)
Adenovirus antigen	6	15.6 ng/ml

TCID₅₀/mL = 50% tissue culture infectious dose

Cross reactivity study

The cross reactivity of the Vstrip Adenovirus Rapid Test was assessed by testing the following organisms and viruses. None of the microorganisms tested in the following table gave a positive result in the Vstrip Adenovirus Rapid Test.

Bacterial Panel

<i>Pseudomonas aeruginosa</i>	<i>Serratia marcescens</i>
<i>Staphylococcus aureus</i>	<i>Staphylococcus epidermidis</i>
<i>Streptococcus Group A</i>	<i>Streptococcus Group B</i>
<i>Streptococcus Group C</i>	<i>Streptococcus Group F</i>
<i>Streptococcus Group G</i>	<i>Streptococcus mutans</i>
<i>Streptococcus pneumoniae</i>	<i>Streptococcus sanguis</i>
<i>Salmonella typhi</i>	<i>Salmonella enteritidis</i>
<i>Vibrio cholera</i>	

Viral Panel

Human respiratory syncytial virus	Human rotavirus
Human parainfluenza virus type 2	Echovirus type 11
Herpes simplex virus type2	

Interfering Substances

Common the-counter medication, anti-diarrhea drug and substances occasionally present in feces were tested with the Vstrip Adenovirus Test in order to verify that they did not interfere with the test.

Interference substances/ Concentration

Aspirin/ 20mg/ml	Metronidazole/ 0.25 mg/ml
Barium sulfate/ 0.25 mg/ml	Voren Supp./ 12.5mg/ml
Bilirubin/ 0.25 mg/ml	PECOLIN/ 5%
Hemoglobin/ 0.25 mg/ml	

PACKAGING

Vstrip Adenovirus Rapid Test.....20 Tests/Kit

ORDERING INFORMATION

Product No. : IGO5020S

REFERENCES

1. David A. Warrell, Timothy M. Cox, and John D. Firth Oxford Textbook of Medicine (5 ed.)
2. Chetana Vaishnavi infections of the Gastrointestinal System
3. Biosafety in Microbiological and Biomedical Laboratories, 5th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (2007).
4. Henretig F.M. MD, King C. MD, Textbook of Pediatric Procedures, Chapter 123 – Obtaining Biologic Specimens Williams and Williams (April 1997).
5. The Clinical Virology Laboratory, Department of Laboratory Medicine at Yale: <http://info.med.yale.edu/labmed/virology/booklet.html>.

SYMBOL LEGEND

Catalog Number	Batch Code	Contains sufficient for < n > tests
In vitro diagnostic medical device	Positive Control	Do not reuse
Consult instructions for use	Use-by date	CE Marking
Temperature Limit	Manufacturer	Authorized Representative in the European Community

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Panion & BF Biotech Inc.

16F., No.3, Park St., Nangang Dist., Taipei City 115, Taiwan

Manufactured by: Panion & BF Biotech Inc. at Xizhi

Manufacturer Address:
6F.-3, No.306, Sec. 1, Datong Rd., Xizhi Dist., New Taipei City 22146 Taiwan

Technical Support: +886-2-2691-9895

Email: xizhi@pbf.com.tw

www.pbf.com.tw

www.vstriptechn.com

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HiSS Diagnostics GmbH, Tullastraße 70 D-79108 Freiburg