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Vstrip Flu A&B Rapid Test.....20Tests/Kit

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


- Murphy B.R. and Webster R.G., 1996, Orthomyxoviruses, pp. 1397–1445. In: Fields Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M.Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 5th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (2007).
- Henretig F.M. MD, King C. MD, Textbook of Pediatric Procedures, Chapter 123 – Obtaining Biologic Specimens Williams and Williams (April 1997).
- The Clinical Virology Laboratory, Department of Laboratory Medicine at Yale: http://info.med.yale.edu/labmed/virology/booklet.html.

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The Vstrip Flu A&B Rapid Test is for the detection of Influenza type A and type B viral antigens in the nasal secretions.

INTENDED USE

Vstrip Flu A&B Rapid Test is a rapid *in vitro* immunochromatographic assay for the qualitative detection of influenza type A and type B virus from the nasopharyngeal swab of symptomatic patients. The influenza A & B virus can be detected by this test in the very early stages of infection when symptoms have just initiated. The test is intended for professional and laboratory use as an aid in the rapid diagnosis of influenza type A and type B viral infections. The test doesn't detect Influenza type C and negative test results should be confirmed by cell culture.

SUMMARY AND EXPLANATION

Influenza viruses are single-strand RNA viruses and include type A, type B and type C. Type A viruses are classically more prevalent and are associated with the most serious influenza epidemics, while type B viruses are generally milder than type A. Type C viruses have not been previously associated with any serious epidemic diseases in history. Both types A and B viruses can prevalent simultaneously, but usually one type is dominant during a given season.¹ Influenza is a highly contagious, acute, viral infection of the upper respiratory tract. Typical symptoms of the influenza infection are an abrupt onset of fever, headache, myalgia, malaise, nonproductive cough, sore throat and rhinitis. The infection usually continues for about a week and most people recover within 1 to 2 weeks. If it occurs to the young, the elderly and people suffering from medical conditions such as lung diseases, diabetes, cancer, kidney or heart problems, influenza poses a serious threat. Under normal circumstances, the infection may lead to severe complications of underlying diseases, pneumonia and death. Vstrip Flu A&B Rapid Test is an immunochromatographic test that using monoclonal antibody-coated latex to detect the presence of Influenza type A and type B viral antigens in the nasal secretions. The test is easy to perform and test results can be visually interpreted in 10 minutes.

PRINCIPLE OF THE TEST

Vstrip Flu A&B Rapid Test is a sandwich solid phase immunochromatographic assay. To perform the test, insert the test dipstick into the extraction buffer. The sample flows through a labeled pad containing anti-influenza A antibodies coupled to blue latex and anti-influenza B antibodies coupled to red latex. If the sample contains influenza A and/or influenza B viral antigens, the antigen will bind to the antibody to form antigen-antibody-latex complexes. These complexes move on the nitrocellulose membrane via capillary action toward the test line region in which influenza specific antibody is spotted. As the complexes reach the test line, they will bind to the antibody on the membrane in the form of a line. A purple control line will always appear in the result window to indicate that the test has been performed correctly and that the test device has functioned properly. If influenza antigen is not present, or is present at very low levels in the specimen, only the control line will be visible. If the control line does not develop, the test is invalid.

MATERIALS AND REAGENTS PROVIDED

- Test Dipsticks (20):
Mouse monoclonal anti-influenza A and anti-influenza B antibodies.
- Extraction Buffer(20):
Tube vials with detergent, protein and salt.
- Nasopharyngeal Swabs (20)
- Influenza A Positive Control Swab (1): Swab is coated with non-infectious recombinant influenza A antigen.
- Influenza B Positive Control Swab (1): Swab is coated with non-infectious recombinant influenza B antigen
- Package Insert (1)

MATERIALS NOT PROVIDED

- Specimen containers
- Timer or watch
- Gloves

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Directions should be read and followed carefully.
- Do not use the kit contents beyond the expiration date printed on the outside of the box.
- Do not interchange or mix different lots of Vstrip Flu A&B Rapid Test.
- Do not insert the test dipstick directly into the sampling area (month, nasal).
- Positive control swab should be mixed in the extraction buffer gently before use.
- Disregard test results beyond specified time (10 min).
- Use appropriate precautions in the collection, handling, storage, and disposal of patient samples and used kit contents.²

- Use of Nitrile or Latex gloves is recommended when handling patient samples.²
- Dispose of containers and used contents in accordance with Local requirements.
- Do not reuse kit components or test devices.
- The test dipsticks must remain sealed in the protective foil pouch until use.
- To obtain accurate results, you must follow the Package Insert.
- Inadequate or inappropriate specimen collection, storage, and transport may yield false test results.
- Seek specific training or guidance if you are not experienced with specimen collection and handling procedures.^{3,4}

STORAGE

- The product should be stored at 15-30°C , away from direct sunlight.
- Kit contents are stable until the expiration date printed on the outer box.
- The test dipsticks must be kept in the sealed pouch until use.
- Do not freeze or overheat the test kit or kit reagents.

SPECIMEN COLLECTION AND HANDLING

Proper specimen collection, storage, and transport are critical to the performance of this test.

Nasal Swab Sample:

It is important to obtain as much secretion as possible. Therefore, to collect a nasal swab sample, insert the sterile swab into the nostril that presents the most secretion under visual inspection. Using gentle rotation, push the swab until resistance is met at the level of the turbinates (less than one inch into the nostril). Rotate the swab a few times against the nasal wall.

Nasopharyngeal Swab Sample:

For optimal test performance, use the swabs supplied in the kit. It is important to obtain as much secretion as possible. Therefore, to collect a nasopharyngeal swab sample, carefully insert the sterile swab into the nostril that presents the most secretions under visual inspection. Keep the swab near the septum floor of the nose while gently pushing the swab into the posterior nasopharynx. Rotate the swab several times.

ASSAY PROCEDURE

All clinical specimens must be at *room temperature* before beginning the assay. **Expiration date:** Check the expiration date on each individual test package or outer box before use. Do not use any test beyond the expiration date on the label.


- Place the swab with sample into the extraction buffer. Roll the swab **ONLY** three (3) times while pressing the head against the bottom and the side of the extraction buffer. Leave the swab in the extraction buffer for 1 minute.
- Roll the swab head against the inside of the extraction buffer as you remove it. Dispose of the used swab in accordance with your biohazard waste disposal protocol.
- Place the test dipstick into the extraction buffer with the arrows on the test dipstick pointing down. Do not handle or move the test dipstick until the test is complete and ready for interpretation.
- Read result at 10 minutes. Some positive results may appear sooner. Do not read the result after 10 minutes.

INTERPRETATION OF RESULTS

Positive result

At 10 minutes, the appearance of ANY shade of a blue and/or red Test Line below the purple Control Line, and the appearance of a purple procedural Control Line indicates a positive result for the presence of influenza A and/or B viral antigen.

- If only a **BLUE** line appears below the Control Line (C), the test result is positive for type A (A).
- If only a **RED** line appears below the Control Line (C), the test result is positive for type B (B).
- If both a **BLUE** and **RED** line appear below the Control Line (C), the test results are positive for both types A (A) and type B (B) .

Negative result

At 10 minutes, the appearance of **ONLY** the purple procedural Control Line indicates influenza A and B antigen were not detected. A negative result indicates that the sample is negative for antigen or the antigen level is below the detection limit.

Invalid

If at 10 minutes, the purple procedural Control Line does not appear, even if a Blue or Red Test Line (A or B) appears, the result is considered invalid. If the test is invalid, a new test should be performed with a new patient sample and a new test dipstick.

QUALITY CONTROL

The Vstrip Flu A&B Rapid Test contains external controls within the kit in order to allow users to confirm the integrity of the test dipstick.

Positive Control

- Place the positive control swab into the extraction buffer.
- Mix the solution by rotating the positive control swab forcefully against the side of the extraction buffer at least 1 minute. Extract as much liquid as possible from the swab, squeezing the sides of the extraction buffer as the swab is withdrawn. Discard the swab with biohazard waste disposal protocol.
- Remove the test dipstick from its protective foil pouch before using it.
- Put the test dipstick into the extraction buffer (the arrows pointing down).
- Read the result within 10 minutes. Do not read the test result later than 10 minutes.

LIMITATIONS

- The contents of this kit are to be used for the qualitative detection of influenza A and B antigen from the Nasal and Nasopharyngeal swab.
- A negative test result may occur if the level of antigen in a sample is below the detection limit of the test.
- Failure to follow the Test Procedure and interpretations of Test Results may adversely affect test performance and/or invalidate the Test Results.
- Test Results must be evaluated in conjunction with other clinical data available to the physician.
- Negative test results do not rule out other potential non-influenza viral infections. Negative results should be confirmed by cell culture.
- Positive test results do not rule out co-infections with other pathogens.
- Positive test results do not identify specific influenza A or B virus subtypes.
- Positive and negative predictive values are highly dependent on prevalence. False negative test results are more likely during peak activity when the prevalence of disease is high. False positive test results are more likely during periods of low influenza activity when prevalence is moderate to low.
- Individuals who received nasally administered influenza vaccine may have positive test results for up to three days after vaccination.
- Monoclonal antibodies may fail to detect, or detect with less sensitivity, influenza viruses that have undergone minor amino acid changes in the target epitope region.
- If differentiation of specific influenza A or B subtypes and strains is needed, additional testing, in consultation with the state or local public health department, is required.
- Background does not affect the interpretation of the results.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (Limit of Detection)

The analytical sensitivity was evaluated using 20 influenza strains; 17 influenza A and 3 influenza B strains. Each strain was serially diluted in Vstrip extraction buffer. Strains were assayed in triplicate using the Vstrip Flu A&B Rapid Test until no positive signal could be observed. Results are summarized below:

The lowest TCID₅₀/mL obtained for Flu A was 2.64 x 10⁵

The lowest TCID₅₀/mL obtained for Flu B was 3.71 x 10⁵

| Viral Strain | Type | Subtype | ATCC No. |
|---------------------------|------|---------|----------|
| Taiwan / 141 / 2002 | A | H1N1 | N/A |
| Taiwan / 2235 / 2009 | A | H1N1 | N/A |
| Taiwan / 2651 / 2010 | A | H1N1 | N/A |
| Taiwan / 12 / 2011 | A | H1N1 | N/A |
| Taiwan / 2030 / 2014 | A | H1N1 | N/A |
| NWS / 1933 | A | H1N1 | VR-219™ |
| PR / 8 / 1934 | A | H1N1 | VR-1469™ |
| FM / 1 / 1947 | A | H1N1 | VR-97™ |
| Taiwan / 3277 / 2011 | A | H3N2 | N/A |
| Taiwan / 2722 / 2012 | A | H3N2 | N/A |
| Taiwan / 2109 / 2013 | A | H3N2 | N/A |
| Taiwan / 2195 / 2014 | A | H3N2 | N/A |
| Taiwan / 534 / 2015 | A | H3N2 | N/A |
| Hong Kong / 8 / 1968 | A | H3N2 | VR-1679™ |
| Taiwan / 1050 / 2014 | B | | N/A |
| Taiwan / 2578 / 2014 | B | | N/A |
| Taiwan / 3859 / 2014 | B | | N/A |
| New Caledonia / 20 / 1999 | A | H1N1 | N/A |
| Anhui/1/2013 | A | H7N9 | N/A |
| Hong Kong / 486 / 1997 | A | H5N1 | N/A |

TCID₅₀/mL = 50% tissue culture infectious dose

Detection Specificity

Vstrip Flu A&B Rapid Test detects the influenza A and/or influenza B virus at >99% specificity.

Cross reactivity study

The Vstrip Flu A&B Rapid Test was evaluated with a total of 15 bacteria, 7 viruses. Bacteria were evaluated at 10⁷ org/ml. Viruses were evaluated at a concentration of at least 10⁵–10⁷ TCID₅₀/mL. None of the microorganisms tested in the following table gave a positive result.

| Bacteria panel | Viral panel |
|-------------------------------|-----------------------------------|
| <i>Pseudomonas aeruginosa</i> | Human respiratory syncytial virus |
| <i>Serratia marcescens</i> | Human rotavirus |
| <i>Staphylococcus aureus</i> | Adenovirus 7 |

| | |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|
| <i>Staphylococcus epidermidis</i> <i>Streptococcus Group A</i> <i>Streptococcus Group B</i> <i>Streptococcus Group C</i> <i>Streptococcus Group F</i> <i>Streptococcus Group G</i> <i>Streptococcus mutans</i> <i>Streptococcus pneumoniae</i> <i>Streptococcus sanguis</i> <i>Salmonella typhi</i> <i>Salmonella enteritidis</i> <i>Vibrio cholerae</i> | Adenovirus 40 human parainfluenza virus type 2 Echovirus type 11 herpes simplex virus type2 |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|

INTERFERING SUBSTANCES

Whole blood, Swinin nasal spray products and common chemicals were evaluated and did not interfere with the Vstrip Flu A&B Rapid Test at the levels tested above:

| Interference substances | |
|------------------------------------------|-------------------------------|
| Whole Blood (4%) | Oxymetazoline HCl (10 mg/mL) |
| Aspirin/ Acetylsalicylic Acid (20 mg/mL) | Phenylephrine HCl (100 mg/mL) |
| Dextromethorphan (10 mg/mL) | Swinin nasal sprays (10%) |
| Diphenhydramine HCl (5 mg/mL) | |

PRECISION STUDIES

The total, within-run, and between-run performance of the Vstrip Flu A&B Rapid Test was evaluated for precision. A panel consisting of two different levels of influenza A antigen (Texas 1/77(H3N2); weak positive and strong positive) and two different levels of influenza B antigen (Hong Kong 5/72; weak positive and strong positive) were repeated five times with a single lot of Vstrip Flu A&B Rapid Test on three different days. >99% accuracy was obtained from all specimens tested.

PACKAGING

Vstrip Flu A&B Rapid Test.....20 Tests/Kit

ORDERING INFORMATION

Product No. : IG01020S

REFERENCES

- Murphy B.R. and Webster R.G., 1996, Orthomyxoviruses, pp. 1397–1445. In: Fields Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M.Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 5th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (2007).
- Henretig F.M. MD, King C. MD, Textbook of Pediatric Procedures, Chapter 123 – Obtaining Biologic Specimens Williams and Williams (April 1997).
- The Clinical Virology Laboratory, Department of Laboratory Medicine at Yale: <http://info.med.yale.edu/labmed/virology/booklet.html>.

SYMBOL LEGEND

| | | |
|-----------------------------------------------|----------------------------------|-------------------------------------------------------------------|
| REF Catalog Number | LOT Batch Code | Contains sufficient for < n > tests |
| IVD In vitro diagnostic medical device | CONTROL+ Positive Control | Do not reuse |
| Consult instructions for use | Use-by date | CE CE Marking |
| Temperature Limit | Manufacturer | EC REP Authorized Representative in the European Community |

07/16

CE **EC REP** HiSS Diagnostics GmbH, Tullastraße 70
D-79108 Freiburg

Panion & BF Biotech Inc.

16F., No.3, Park St., Nangang Dist., Taipei City 115, Taiwan

Manufactured by: Panion & BF Biotech Inc. at Xizhi

Manufacturer Address:

6F.-3, No.306, Sec. 1, Datong Rd., Xizhi Dist., New Taipei City

22146 Taiwan

Technical Support: +886-2-2691-9895

Email: xizhi@pbf.com.tw

www.pbf.com.tw

www.vstriptechnology.com