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Vstrip Rotavirus Rapid Test for the Detection of Rotavirus Antigen in Human Stool Specimen

INTENDED USE

Vstrip Rotavirus Rapid Test is an in vitro qualitative immunochromatographic assay for the detection of Rotavirus antigens in human stool specimen. Test results are intended to aid in the early diagnosis of Rotavirus infection or early confirmation of infection. This test is for professional use only. If the check result is uncertain, and patients with significant clinical symptoms, a further determination should be carried out with other clinical diagnostics.

SUMMARY AND EXPLANATION

Rotavirus is the major cause of children with acute gastroenteritis, especially in children 6 to 24 months in age^{1,3}. The average incubation period for the Rotavirus infection is 3 days. Initial symptoms of Rotavirus infection are fever and vomiting, followed by frequent watery diarrhea, and diarrhea usually lasts 3-9 days. Severe dehydration can cause death. 2 days before symptoms appear to 10 days after finding a large number of Rotavirus in feces. The gold standard for detecting Human Rotavirus depends on electron microscopy since existing tissue culture methods are unreliable³. However, newly introduced enzyme immunoassays and latex agglutination assays with increased sensitivities and specificities are now the methods of choice^{2,4}. Vstrip Rotavirus Rapid Test offers a simple, rapid method for detecting Rotavirus antigen in patient stool.

PRINCIPLE OF THE TEST

Vstrip Rotavirus Rapid Test is a qualitative immunochromatographic assay for the determination of Rotavirus in stool samples. The membrane is pre-coated with monoclonal antibodies, on the test band region, against viral antigens. During testing, the sample is allowed to react with the colored conjugate (anti-Rotavirus monoclonal antibodies) which was pre-dried on the test strip. The mixture then moves upward on the membrane by capillary action.

An extract is first prepared by suspension of the specimen in the provided diluent buffer solution. If the sample extract contains Rotavirus antigens, these will form an antigen-antibody complex with the colored particles and a visible Test band will appear as a Red line. If Rotavirus antigen is not present, or is present at very low level in the specimen, only the control line will be visible. The C band (Blue color) should always appear after a sample is applied, indicating a valid result. By this means, the device can accurately indicate the presence of Rotavirus in 10 minutes.

MATERIALS AND REAGENTS PROVIDED

- Test dipsticks (20)
Each strip contains Rotavirus specific monoclonal antibodies.
- Stool collection tubes with diluent (20)
The buffer contains 0.08% sodium azide.
- Rotavirus Positive control swab (1)
- Package Insert (1)

MATERIALS NOT PROVIDED

- Specimen collection container
- Timer
- Disposable gloves

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use.
- Read the entire procedure carefully prior to testing.
- Do not open the test pouches until ready to perform the assay.
- Use of Nitrile or Latex gloves is recommended when handling patient samples⁵.
- Stool must be mixed thoroughly to insure a representative sample prior to pipetting. Do not use stools that have dried out.
- Patient specimens should be mixed in the stool collection tubes with diluent gently and thoroughly before use.
- Disregard test results beyond specified time (10 min).
- Do not insert the test strip directly into patient stool specimens.
- Use only one sample diluent vial. Discard after use. Do not attempt to reuse.
- Do not use the test if the foil pouch is damaged.
- Do not reuse kit components or test devices.
- Do not mouth pipette samples or reagents.
- Patient specimens may contain infectious agents and should be handled and disposed of as potential biohazards.
- Do not interchange or mix different lots of Vstrip Rotavirus Rapid Test.
- Do not use the kit components beyond the expiration date, printed on the outside of the box.
- Some reagents in this kit contain sodium azide, which is a skin irritant. Avoid skin contact with reagents.
- Dispose all used materials in the appropriate container. Treat as potentially biohazard⁵.
- To obtain accurate results, you must follow the Package Insert.

SHELF LIFE AND STORAGE INSTRUCTION

The expiration date is indicated on the package label. Store kit reagents at 15-30°C, away from direct sunlight.

SPECIMEN COLLECTION AND STORAGE

Proper specimen collection, storage, and transport are critical to the performance of this test.

For the best results, specimens should be collected after onset of symptoms. Do not leave specimens at room temperature for prolonged periods. Specimens containing high levels of blood may fail to flow in the Vstrip Rotavirus device, resulting in an invalid / false positive test result. Testing of an additional specimen is recommended under such circumstances. Stool specimens may be collected by the method routinely utilized by the laboratory, providing no dilution of the specimen occurs.

Stool specimens should be collected in clean containers that do not contain media, preservatives, animal serum, or detergents as any of these additives may interfere with the Rotavirus antigen test or obtained from a diaper. Specimens may be stored at 2~8°C for 2 days without interfering with the test performance. For long-term storage of specimens, -20°C or colder is recommended. Repeated freezing and thawing of specimens is not recommended and may cause erroneous results. Do not store specimens in self-defrosting freezers.

SPECIMEN PREPARATION

Stool specimen:

- Unscrew the sample diluent vial cap, use the spoon attached to the cap to transfer small piece of stool (2-5 mm in diameter; approximately 30-50 mg) into the vial containing sample diluent buffer.
- Close the vial and tighten securely. Shake the stool collection tube in order to assure good sample dispersion.
- Label the sample diluent vial with a patient or other control information.

ASSAY PROCEDURE

- Remove a test strip from its sealed bag.
- Unscrew the sample diluent vial cap then place the test strip into the stool collection tube (the arrows pointing down).
- Read the result at 10 minutes. Some strong positive results may appear sooner. Do not read result after 10 minutes.

INTERPRETATION OF RESULTS

Positive Result:

Visually detectable in addition to the blue colored line on the control region, a distinct red colored line also appears in test region. A positive result indicates the presence of Rotavirus antigen.

Negative Result:

Visually detectable no line appears in the test region. A distinct blue line shows on the control region. A negative result indicates that Rotavirus antigen is absent or below the level of detection.

Invalid:

No visually detectable the blue colored line in the control region, with or without a visually detectable red colored line on the test region. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

Notes on the Interpretation of Results


The intensity of the red colored band in the result will vary depending on the concentration of antigens in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

The Positive Controls should be assayed once upon receipt of the kit. Internal procedural controls are included in the test. A blue line appearing in the control region is an internal control. It confirms sufficient specimen volume and correct procedural technique. The positive controls are intended to monitor for reagent failure, but will not ensure precision at the analytical assay cutoff.

QUALITY CONTROL TESTING PROCEDURE

Positive Control:

- Unscrew the sample diluent vial cap and immediately introduce the positive control swab into the stool collection tube.
- Mix the solution by rotating the positive control swab forcefully against the side of the stool collection tube at least 1 minute. Extract as much liquid as possible from the swab, squeezing the sides of the stool collection tube as the swab is withdrawn.
- Remove the test strip from its sealed bag before using it.
- The test strip put into the stool collection tube (the arrows pointing down).
- Read the result within 10 minutes. Do not read the test result later than 10 minutes.
- Result: The Positive Control should yield visually detectable two lines (Blue: control line/Red: test line).

LIMITATIONS OF THE PROCEDURE

- Vstrip Rotavirus Rapid Test is highly sensitive and specific for Rotavirus antigen. The monoclonal antibody in this test will detect human Rotaviruses, but cannot be used to differentiate types. As with all in vitro diagnostic procedures, test results should be interpreted by a physician in conjunction with other clinical information.
- After one week of infection, the number of viruses in feces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
- An excess of stool sample could cause wrong results (brown bands appear).
- A negative result does not exclude the possibility of Rotavirus infection in the patient. False negative results may occur due to low concentration levels of the Rotavirus antigen below the sensitivity level of the test, improper or inadequate sampling, or improper handling of the specimen. Test results should be interpreted by a physician
- Co-infection with bacterial pathogens is possible. Therefore, bacteriological tests should be performed in parallel with this test to rule out bacteriological etiology.

EXPECTED VALUES

Vstrip Rotavirus Rapid Test detects the presence of Rotavirus antigen in stool. Expected values for a given population should be determined for each laboratory. The prevalence of Rotavirus infection will vary based on many factors such as age, geographic location, method of sample collection, sample handling and transportation, and the general health environment of the patient population under study. Healthy individuals tested should be negative for Rotavirus. Some infected individuals may show symptoms or only minor symptoms, and these patients may test negative.

PERFORMANCE CHARACTERISTICS

The evaluation was conducted comparing the results obtained using Vstrip Rotavirus Rapid Test to another commercially available Rotavirus Rapid Test .

Sensitivity

The detection of Rotavirus showed >99% of concordance in sensitivity.

Specificity

The detection of Rotavirus showed > 99% of concordance in specificity.

Reproducibility

Three different lots of the Vstrip Rotavirus Rapid Test was tested in multiple replicates using negative, low positive, medium positive and high positive samples. Reproducibility was tested within each lot, between lots, and on different days. In each case, all tests yielded >99% reproducibility.

Limit of Detection

To determine the limits of detection of the Vstrip Rotavirus Rapid Test. The study used human Rotavirus (HRV) as following table which were obtained from ATCC. The limit of detection of the test is a concentration of 1.58x10³ TCID₅₀/mL.

Rotavirus strain	Limit of detection (TCID ₅₀ /mL)
Hu/Australia/102510/77/L	8.89X10 ²
HRV 89-12C2	1.58X10 ³

Cross Reactivity

The cross reactivity of the Vstrip Rotavirus Rapid Test was assessed by testing the following organisms and viruses. None of the microorganisms tested in the following table gave a positive result in the Vstrip Rotavirus Rapid Test.

Microorganisms

<i>Pseudomonas aeruginosa</i>	<i>Serratia marcescens</i>
<i>Staphylococcus aureus</i>	<i>Staphylococcus epidermidis</i>
<i>Streptococcus Group A</i>	<i>Streptococcus Group B</i>
<i>Streptococcus Group C</i>	<i>Streptococcus Group F</i>
<i>Streptococcus Group G</i>	<i>Streptococcus mutans</i>
<i>Streptococcus pneumoniae</i>	<i>Streptococcus sanguis</i>
<i>Salmonella typhi</i>	<i>Salmonella enteritidis</i>
<i>Vibrio cholera</i>	

Viruses

Adenovirus type 7 Adenovirus type 41
Echovirus Type 11

Interfering Substances

Common the-counter medication, anti-diarrhea drug and substances occasionally present in feces were tested with the Vstrip Rotavirus Test in order to verify that they did not interfere with the test. The following substances showed no interference:

Interference substances/ Concentration

Aspirin/ 20mg/ml	Metronidazole/ 0.25 mg/ml
Barium sulfate/ 0.25 mg/ml	Voren Supp./ 12.5mg/ml
Bilirubin/ 0.25 mg/ml	PECOLIN/ 5%
Hemoglobin/ 0.25 mg/ml	

PACKAGING

Vstrip Rotavirus.....20 Tests/Kit

ORDERING INFORMATION

Product No. : IG04020S

REFERENCE

- Burke, B. and U. Desselberger. 1996. Minireview: Rotavirus Pathogenicity. *Virology* 218: 299-305.
- Christensen, M.L. 1989. Human viral gastroenteritis. *Clin. Microbiol. Rev.* 2: 51-89.
- Truant, A.L. and T. Chonmaitree. 1982. Incidence of Rotavirus infection in different age groups of pediatric patients with gastroenteritis. *J. Clin. Microbiol.* 16: 568-569.
- Mathewson, J., Winsor, D. Jr., Dupont, H.L., and Secor, S. L. (1989). Evaluation of assay systems for the detection of Rotavirus in stool specimens. *Diag. Microbiol. Infect. Diseases.* 12:139-141.
- Biosafety in Microbiological and Biomedical Laboratories, 5th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (2007).

SYMBOL LEGEND

Catalog Number	Batch Code	Contains sufficient for < n > tests
In vitro diagnostic medical device	Positive Control	Do not reuse
Consult instructions for use	Use-by date	CE Marking
Temperature Limit	Manufacturer	Authorized Representative in the European Community

07/16

HiSS Diagnostics GmbH, TullastraÙe 70
D-79108 Freiburg

Panion & BF Biotech Inc.

16F., No.3, Park St., Nangang Dist., Taipei City 115, Taiwan

Manufactured by: Panion & BF Biotech Inc. at Xizhi

Manufacturer Address:

6F.-3, No.306, Sec. 1, Datong Rd., Xizhi Dist., New Taipei City 22146 Taiwan

Technical Support: +886-2-2691-9895

Email: xizhi@pbf.com.tw

www.pbf.com.tw

www.vstriptechnology.com